Detektor bublin (Air Bubble Detector)

Typ SONOCHECK **ABD05**

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Obsah

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Událost / Operating mode	Priorita	Signál na ABD-OUT
Kapalina	Low	L
Vzduch / Bublina	Medium	H
Interní chyba (self-test)	High	H
Závažná porucha (watchdog)	Highest	H
Start procedure after power on or reset (approx. 1 s)	--	H Note: After 'Init Interval' the 'Initial Test' will be performed. This includes test routines for checking output ABD-OUT. Level of output will change several times.
Boot mode	--	UART Tx Note: Serial data are transmitted on demand. In inactive state = H

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The start pattern and the CRC guarantee that the signal sequence is interpreted accurately.

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
Událost / Operating mode	Priorita	Data na ABD-OUT (Tx portu UART)
Kapalina	Low	0 ... 239, úměrně k velikosti bubliny
Vzduch / Bublina	Medium	240
Interní chyba (self-test)	High	241
Serious fault (watchdog)	Highest	Static High
Start procedure after power on (approx. 1 s)	--	High Note: After 'Init Interval' the 'Initial Test' will be performed. This includes test routines for checking output ABD-OUT. Level of output will change several times.
Boot mód	--	UART Tx Note: Serial data are transmitted on demand. In inactive state = High

Tabulka 4: Výstup v Interface módu 1 / 2

To check the sensor functionality, the serial command `CmdBubbleTest` can be used.

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In this mode the sensor works in dialog mode. Information and data are only sent on request by the machine. The interface mode 3 is not suitable for bubble detection in real time.

5.5 Interface mód 4 / 5

Výstup je nastaven jako spínací, vstup je nastaven jako testovací

'Bezpečný stav' snímače Interface Mode 4: Output ABD-OUT = High
Interface Mode 5: Output ABD-OUT = Low

The assignment to events can be configured by parameters. Default settings are:


Událost / Operating mode	Priorita	Signál na ABD-OUT	
		Interface Mode 4	Interface Mode 5
Liquid	Low	Low	High
Air / Bubble	Medium	High	Low
Internal error (self-test)	High	High	Low
Serious fault (watchdog)	Highest	High	High Note: If the 'Bubble Test' is applied, this state could be reliably recognized as fault by the machine.
Start Mode	--	High Note: Note: After 'Init Interval' the 'Initial Test' will be performed. This includes test routines for checking output ABD-OUT. Level of output will change several times.	
Boot mode	--	UART Tx Note: Serial data are transmitted on demand. In inactive state = High	

Tabulka 4 Výstup v Interface módu 4 / 5

⚠ WARNING To reach a high level of safety it is essential, that the output is checked during operation with filled tubing. This can be achieved by periodically performing the 'Bubble Test'.

Operating Mode	Specification of Input ABD-IN
Normal Mode	The input is configured as test input to release a 'Bubble Test'. Interface Mode 4: Low-active pulse Interface Mode 5: High-active pulse For the duration of active pulse the sensor simulates a bubble. The output will react accordingly, and the machine can check the outputs. The min. pulse width is 300 µs.
Start Mode	The input is configured after power on or after restart during 'Init Interval' as serial input Rx of the UART. The 'Init Interval' takes approx. 1 s. During this period the sensor reacts to all commands. In such a way the boot mode can be started in every sensor configuration.
Boot Mode	UART Rx Note: Serial data are transmitted on demand. In inactive state = High

Table 5: Vstup v Interface módu 4 / 5

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'Bezpečný stav' snímače The 'Safe State' is defined as one of the following states:
 ↳ ABD-OUT = static High or static Low
 ↳ The modulation of PWM output is set to the default values for bubble (80 %) or fault (90 %).
Note: The values are configurable.

The output signal can be configured by parameters. The default settings are:


Událost / Operating mode	Priority	Signál na ABD-OUT
Liquid / small bubbles	Low	PWM signal: 20 ... < 80 %, proportional to the bubble size
Air / Bubble	Medium	80 %
Internal error (self-test)	High	90 %
Serious fault (watchdog)	Highest	Static High
Start Mode	--	High Note: After 'Init Interval' the 'Initial Test' will be performed. This includes test routines for checking output ABD-OUT. Level on output will change several times.
Boot mode	--	UART Tx Note: Serial data are transmitted on demand. In inactive state = High

Tabulka 6: Vstup v Interface módu 6 / 7

Operating Mode	Specification of Input ABD-IN
Normal Mode	The input is configured as test input to release a 'Bubble Test'. Interface Mode 6: active low pulse Interface Mode 7: active high pulse For the duration of the active pulse the sensor simulates a bubble. The output will react accordingly and the machine can check the outputs. The min. pulse width is 300 µs.
Start Mode	The input is configured after power on or after restart during 'Init Interval' as serial input Rx of the UART. The 'Init Interval' takes approx. 1 s. During this period the sensor reacts to all commands. In such a way the boot mode can be started in every sensor configuration.
Boot Mode	UART Rx Note: Serial data are transmitted on demand. In inactive state = High

Tabulka 7: Vstup v Interface módu Mode 6 / 7

To reach a high level of safety it is essential, that during operation the output is checked with a filled tubing. This can be achieved by periodically performing the 'Bubble Test'.

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'Bezpečný stav' snímače The 'Safe State' is defined as one of the following states:

- ↳ ABD-OUT = static High or static Low
- ↳ The modulation of PWM output is set to the default values for bubble (80 %) or fault (90 %).

Note: The values are configurable.

The output signal can be configured by parameters. The default settings are:

Událost / Operating mode	Priorita	Signál na ABD-OUT
Liquid	Low	PWM signal: 20 ... < 80 %, represents bubble size
Air / Bubble	Medium	80 %
Internal error (self-test)	High	90 %
Serious fault (watchdog)	Highest	Static High
Start Mode	--	High Note: After 'Init Interval' the 'Initial Test' will be performed. This includes test routines for checking output ABD-OUT. Level on output will change several times.
Boot mode	--	UART Tx Note: Serial data are transmitted on demand. In inactive state = High

Tabulka 8: Výstup v Interface módu 8

The Input is configured as serial input of the UART. The sensor is able to react to serial commands.

To reach a high level of safety it is essential, that the output is checked during operation with filled tubing. This can be achieved by periodically sending the command CmdBubbleTest.

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For the use in medical devices (MDEV) special regularities apply. The following information and instructions partly refer to the previous chapters. Please make sure, that you have read and understood the description of the interface and the interface modes.

⚠ WARNING The sensor offers different features to ensure a high level of patient safety:

- ↳ Fail-safe architecture, incl. self-test procedures and 'Initial Test'
- ↳ Secured data transmission via serial interface or PWM
- ↳ Switching output and periodical 'Bubble Test'.

It is in the responsibility of the manufacturer of the medical device (MDEV) to establish a safe and efficient interface between sensor and MDEV.

The sensor provides safety related information to the MDEV. As the sensor has no direct access to controlling components, like e.g. pumps or valves, it is the responsibility of the MDEV to read all safety related information the sensor provides and to react in an appropriate manner.

The SONOCHECK ABD05 together with the MDEV and other components of the MDEV (like e.g. pumps or valves) form a system for the protection of the patient. An alarm handling should be implemented in the medical device.

The sensor provides different options to achieve a high degree of safety. Via software different input/output configurations for the interface between sensor and machine can be selected. The following chapters describe these different options

6.1 'Fail-Safe architektura

The sensor is connected to the controlling and protective unit of the medical device (MDEV). It is classified under the Medical Device Directive 93/42/EEC as a class IIb product. The sensor design includes a 'Fail-Safe-Architecture'.

Please be aware, that a 'fail safe' sensor might fail, but this will always happen in a safe and secure way.

The sensor performs a 'Cyclical Self-Test' during operation. Furthermore, a black box test is implemented, the 'Initial Test'. This 'Initial Test' is executed after every restart or power-on.

The sensor supports different interfaces to the MDEV:

- Switching output and digital input ABD-IN ('Bubble Test'), to check the proper function of the output. The MDEV should periodically send a test signal. Thus the correct function of the switching output can be verified by the MDEV.
- A Serial interface for a real-time transmission of data, including the bubble size and alarm state to the MDEV, and for controlling the sensor.
- Pulse-Width-Modulated output signal, with coded information about the bubble size.

6.2 'Initial Test' – vstupní boot mód

After Power-on or reset the sensor starts with the 'Init Interval'. It takes approx. 1 s. Within this period the sensors reacts to serial commands. In such a way the sensor can always be switched to the boot mode, independent of the interface configuration.

The 'Initial Test' follows the 'Init Interval'. This 'Initial Test' includes the test of the watchdog and of the output ABD-OUT. Sleeping failures or other second failures will be recognized by the sensor during this test.

6.3 'Bubble Test'

The test is implemented to check the correct function of the output. The test has to be initiated by the medical device. The procedure depends on the interface mode:

Mód 4 / 6: For the duration of the active low signal on input ABD-IN, the sensor simulates a bubble.

Mód 5 / 7: For the duration of the active high signal on input ABD-IN, the sensor simulates a bubble.

Mód 0 / 1 / 2 / 8: For the duration of a defined cycle count the sensor simulates a bubble. The cycle count is a parameter of the sensor. It can be modified by a serial command. Each cycle equals 200 μ s.

The sensor will react to the simulated bubble and set the output ABD-OUT accordingly. The output should be checked by the MDEV:

Mód 0 / 4: The switching output will be set to High for the period of the simulated bubble

Mód 5: The switching output will be set to Low for the period of the simulated bubble

Mód 1 / 2: The transmitted serial data will indicate the simulated bubble

Mód 6 / 7 / 8: The PWM signal will indicate the simulated bubble

For the bubble simulation the amplitude of the transmitted pulse is reduced by approximately – 8 dB. This reduced signal is processed by the sensor in the same way as a real bubble would be processed. That means, the sensor does not differentiate between a real bubble and the test.

To ensure, that the sensor is working properly, the medical device should trigger this bubble test periodically, when the sensor is in the state 'Liquid' – which means the sensor is working properly and has not detected real bubbles. In all other states the test has no effect, since the sensor is already in a 'Safe State'.

6.4 Čas odezvy

The internal measuring cycle is 200 μ s. Each attenuation of the received acoustic signal is evaluated. To avoid false alarms smoothing functions are implemented.

Switching Output

In each cycle of 200 μ s the sensor evaluates the maximum attenuation of the received signal. In the default configuration of the sensor, the response time to bubbles is approx. 1 ... max. 2 ms. With special settings (like e.g. no smoothing function) the response time can be reduced down to 0.2 ms.

Serial Output / PWM Output

The output interval via the serial interface or PWM period is 5 cycles (=1 ms). This results in a response time from the detection of a bubble of 1 to max. 2 ms. A further delay and a dwell time can be configured by parameters.

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WARNING

There is no optical or acoustical alarm. The sensor provides LEDs to signalize the state of bubble detection only for information. This LED should not be used as a part of an alarm messaging system.

The LEDs are freely configurable. They can be assigned to internal events, like bubble, fault or liquid, or they can be remote controlled by the MDEV in Interface Mode 0 / 1 / 2 and 8. The LEDs shall support the user during the treatment only for information, it should not be considered as an optical alarm.

The LED's can be set remotely by the medical device or automatically by the sensor itself.

6.6 Watchdog

A watchdog is implemented in the sensor. The watchdog has direct control of the ABD-OUT output signal. When the self-test detects a failure of the software or the hardware, the watchdog sets the output to High, what is the 'Safe State'. If the output is configured as serial output or PWM output, the medical device can recognize this failure, as the continuously sent data or information on the respective output are missing.

Additionally the watchdog serves as voltage monitor for the operating voltage of 3.3 V, supplied to the microcontroller. In the case of overvoltage, the watchdog sets the output line ABD-OUT to High, what is the 'Safe State'.

The full functionality of the watchdog, including the voltage monitor is tested in the 'Initial Test' after power on or restart.

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Item	Meaning
Architecture	Single channel architecture / Fail-safe
Software	Software developed acc. DIN EN 62304:2006, classified as "C" (Module Main).
Security	The sensor is developed according to standards and directives for medical devices: IEC 60601-1:2005.
'Initial Test' Self-test	The sensor performs an 'Initial Test' each time after power on or restart. The sensor has to be restarted by the medical device at least once a day, to ensure regular execution of the 'Initial Test'. During normal operation, the self-test is performed at least one time in each 0.3 s time frame.
EMC	The sensor was designed to fulfill the requirements according to IEC 60601-1-2: 2007 (3rd edition). This was verified at SONOTEC by realistic tests. The test should be performed by the manufacturer of the medical device in mounted state.
Electrical safety	Classified as applied part "CF" in combination with machine and tubing. Insulation between the outside of the housing and electrical components: > 500 V.
Acoustical safety	The sensor fulfills the requirements according to EN 61157:2008, and is suitable for the use with human blood.
Restriction of hazardous Substances	The sensor fulfills the requirements according to the DIRECTIVE 2011/65/EU 8th of June 2011. The sensor is RoHS compliant, with permissible exemptions under EU-RoHS III 7cl and IV 15
Traceability	Traceability is ensured for all relevant parts of the sensor.

Tabulka 9: Safety / directives SONOCHECK ABD05

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The bubble sensitivity of the sensor depends on the diameter of the tubing and on the mounting position of the sensor. Bubbles with a diameter greater than 30 % ... 50 % of the inner diameter of the tubing will be detected.

⚠ WARNING

The mounting position influences the bubble sensitivity. Avoid sharp bends of the tubing within 5 cm before and after the sensor.

**POZNÁMKA!**


For each application the sensor has to be adjusted to the applied tubing carefully! It is recommended to evaluate the bubble sensitivity for the application with the help of the ABD Monitor!

When the sensor is installed as recommended, bubbles with the following size can be detected:

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Tabulka 10: Bubble sensitivity depending on sensor type

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

Bubbles always move along the highest path in the tube.


Bubbles pass the sound path (symbolized by the blue wavy lines) in center of the measuring channel

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NEDOPORUČUJEME!

Bubbles do not pass the sound path (symbolized by the blue wavy lines) in the center of the measuring channel

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Snímač	Upevnění
Rectangular	Using recessed M4 threaded holes on the rear of the sensor
Circular	Using an appropriate clamping fixture

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POZNÁMKA!

During operation, the tubing must be securely held in the channel of the sensor. If this is not possible, an additional fixture should be installed.

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⚠ WARNING

There is no protection against reverse polarity implemented. In the case of reverse polarity the internal circuit is shorted and a high current occurs. Protect the digital input and output against overvoltage! There is no fuse in the sensor. The maximum current must be limited to 200 mA! The power consumption should be limited in device.

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8.1 Operating information

⚠ WARNING

Incorrect installation and use of the ABD05 sensor and its components can constitute a hazard to the user:

- ↳ The sensor and the associated tubing must be clean and intact during the operation.
- ↳ The sensor channel should be kept free of liquids. Otherwise, particularly with small channel widths, an acoustic short circuit of the measuring cell can occur. This might severely impair the functionality of the sensor.
- ↳ It is not permitted to use tools, especially sharp tools (e. g. screw driver or scissors) to press the tubing into the channel. It is only allowed to press the tubing by hands.

Neither coupling gel nor any other liquid is allowed in the channel of the sensor for coupling, as it might cause an acoustic short circuit, which would strongly disturb the functionality of the sensor.

**POZNÁMKA!**

Continuous use and significant temperature fluctuations change the flexibility of the tubing, the propagation of the ultrasonic signal and therefore the acoustic conditions. The sensor automatically adjusts itself to these changes.

8.2 'Initial Test' po zapnutí nebo po restartu

After power on or restart the sensor is ready for measuring within approximately 1.5 seconds

**POZNÁMKA!**

The sensor starts with 'Init Interval', which takes approximately 1 second. Then the 'Initial Test' is performed, which takes max. 0.5 seconds. After that the normal operation mode is started.

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⚠ WARNING

Incorrect cleaning of the ABD05 sensor can constitute a hazard to the user.

Cleaning is prohibited:

- ↳ in a steam sterilizer, with hot steam in general or by ethylene oxide, gamma- or beta-sterilization
- ↳ with white spirit, acetone or acetone-based solutions
- ↳ by immersion in solvents or other liquids.

The sensor should be cleaned using standard commercial cleaning agents. Standard commercial spray disinfectants should be used for disinfection.

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The error handling can be adapted to individual requirements.

The sensor has self-test routines and sets defined output signals if an error is detected. See chapter 5 for details.

The LED in the channel lights up blue, if an error occurs.

The SONOCHECK ABD05 is maintenance-free.